Anästhesie-Maske

REF

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Quantity:
1 pce

PRODUCT DESCRIPTION

The anaesthesia mask is a component of the mechanical breathing device used in the treatment of patients with respiratory failure.

INTENDED USERS

Nursing staff, medical professionals

INDICATIONS

The product is used as a medical unit for one time use connecting the breathing tube to deliver the anesthetic gas when the patient inhales the anesthetic gas.

CONTRAINDICATIONS

• The product should be contraindicated in patients allergic to the product material.

APPLICATION DESCRIPTION

- To inflate mask cushion, squeeze valve „O“ ring and inflate with syringe.
- In an emergency, squeeze valve „O“ ring and inflate with mouth.
- To deflate mask cushion, squeeze valve „O“ ring and allow cushion to deflate.

Placement of Mask:

- Place mask cushion on the patient's chin first.
- Roll cushion up the face with gentle, but firm pressure.
- Refrain from pressing rigid parts of the mask against the face, nose, or eyes.
- Gentle but firm pressure will provide an effective mask seal.
- If mask contains a detachable prong ring, it can be used in conjunction with a headstrap or removed when not needed.

WARNING NOTICES

- Be sure all connections are secure.
- Do not use in the presence of flammable anesthetics. Do not rinse, soak or wash.
- Excessive pressure applied to face mask can lead to facial and/or optic nerve damage.

LOT NUMBER

The LOT number is printed on the packaging.

EXPIRY DATE

The expiry date is imprinted on the product packaging.

STORAGE CONDITIONS

Do not store in direct sunlight, expose to extreme temperatures or humidity. Store in dry conditions, in a non-corrosive atmosphere and in a well-ventilated room.

HANDLING CONDITIONS

In case of damaged packaging or inappropriate storage the product SHOULD NOT be used and has to be destroyed.

STERILE CONDITION

The product is sterilised by ethylenoxide. It is delivered in a single sterile barrier system (○) with a protective outer packaging (○).

SINGLE USE NOTICE RESP. NOTICE ON REUSE

The product is a disposable product. Do not reuse.

REPROCESSING NOTICE

The product SHOULD NOT be reprocessed.

DISPOSAL

The product can be disposed with normal clinical waste or according to local regulations. Special guidelines need to be taken into consideration..

NOTICE ON REPORTING REQUIREMENTS

Please inform the manufacturer, the distributor and the importer in all cases of serious incidents in connection with the product.

 Hangzhou Shanyou Medical Equipment Co., Ltd.
No. 138, Louta Development Zone, Guancun Village,
Louta Town, Xiaoshan District, Hangzhou, Zhejiang
311266, China

 Prolinx GmbH
Brehmstr. 56 | 40239 Düsseldorf, Germany
Tel: +49-2131-4051968-0 | E-mail: info@eulinx.eu

 Importeur
Novo Klinik-Service GmbH
Technologiepark West - Haus 15
Zum Frenser Feld 1
50127 Bergheim | Germany
Tel. +49 2271 99404-0
Fax +49 2271 99404-100
info@its-novo.de | its-novo.de